


Automatic Characterization of Focal Adhesions in TIRF Microscopy Images

Matthew E. Berginski


University of North Carolina at Chapel Hill

7/29/2009 Hahn Lab Meeting



Focal What?


- Points of contact between the substrate and cells
- Consist of dozens of dynamically recruited proteins
- Important in understanding how cells navigate and sample the environment



E. Zamir and B. Geiger. JCS, 114(20):3583-3590, 2001.

Focal What?


- ▶ Points of contact between the substrate and cells
- ▶ Consist of dozens of dynamically recruited proteins
- ▶ Important in understanding how cells navigate and sample the environment



E. Zamir and B. Geiger. JCS, 114(20):3583-3590, 2001.

Paxillin Info/Example Movie

- ▶ Paxillin - adhesion scaffolding protein
- ▶ Recruited to adhesions early, remains associated
- ▶ Data consists of time lapse images of GFP-labeled Paxillin in 3T3 fibroblasts



Outline

Introduction

Finding/Tracking Adhesions

- Identifying Focal Adhesions


- Visualizing the Results

Adhesion Properties




Conclusions

Identification Overview

- ▶ Prior Methods
- ▶ Multistage process
 - ▶ Finding Adhesions
 - ▶ Tracking Adhesions
 - ▶ Visualization



Prior Methods




D. Webb, et al. NCB, 6(2):154-161, 2004.

Determining Adhesion Pixels

- ▶ Technique adapted from E. Zamir, et al. JCS, 112(11):1655-1669, 1999.
- ▶ High-pass image filter applied
- ▶ Threshold used to select pixels




Threshold Variation




Assigning Pixels to Adhesions


- ▶ Technique adapted from E. Zamir, et al. JCS, 112(11):1655-1669, 1999.
- ▶ High-passed pixel values are sorted, then added singly
- ▶ Each added pixel, joins a current adhesion or starts a new adhesion
- ▶ Merging limited to small adhesions





Visualizing Single Adhesions



Tracking Example Movie - Unique Colors





Tracking Example Movie - Birth Time Colors





Properties Overview

- ▶ General Adhesion Lifecycle properties:
 - ▶ Adhesion birth and death locations
 - ▶ Adhesion lifetimes
- ▶ Assembly and Disassembly Rates


General Dynamics Properties





General Dynamics Properties





Single Adhesion Assembly and Disassembly Rates




Single Adhesion Assembly and Disassembly Rates



Rate Fit Qualities




Overall Rate Summary



Conclusions/Future Work

- ▶ System for extracting properties for labeled adhesions:
 - ▶ Finding
 - ▶ Tracking
 - ▶ Visualization
 - ▶ Properties
- ▶ Future Work:
 - ▶ integration of biosensor data



Questions?